



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre

77 Heyford Park

Upper Heyford

Oxfordshire

OX25 5HD

UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for diagnostic or therapeutic use.**

## SB-PDL-100 - Poly D lysine coated glass coverslips (pack of 100)

Size: Pack of 100

### High quality Poly-D-lysine coated glass COVERSLIPS for Immunofluorescence

#### DESCRIPTION

Glass coverslips are coated with poly D lysine on both sides for cell culture and microscopy.

Packing size and dimensions: 100 coverslips, Ø 15 mm x 0.17 (#1.5)

#### SHIPPING AND STORAGE

- Can be shipped at room temperature.
- Can be stored at room temperature for 3 months.
- Store at 4°C for 6 months and at -20°C for long term storage.

Expiration: 1 year at -20°C

#### INTENDED USE

This product is intended for research use only.

#### HANDLING

1. Take a package from freezer and spray with 70% ethanol in tissue culture hood and let it dry completely.
2. Handle the coverslips with sterile sharp forceps while using.
3. Store unused coverslips at room temperature in tissue culture hood and use them within a week.

(If longer storage is desired, store at 4°C).

#### Removing Coverslip from Well after Culture

1. Create a lifting tool by bending a needle 0.5-1 mm from tip.
2. Use one hand with the lifting tool to lift coverslip from its edge and use other hand with forceps to hold and remove coverslip from well.

#### Coating Compatibility

Different cells require different coating substrates. If this type of coating did not produce satisfactory results, try a different coating substrate.