



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

SB-SHIFIX50PLUS - SB-SHIFIX50PLUS Coverslips

Size: Pack of 50


Shi-fix™ coverslips

Do you wish immunofluorescence of suspension cells was as easy as with adherent cells? Now it is! Shi-fix™ (manufactured by Shikhar Biotech) cover-slips allow suspension cells to be layered on them or directly grown as a monolayer. Easy, hassle-free protocol and no spin columns needed for fixing suspension cells to slides. Just add your cells to the coverslips, wait for 30 mins, wash unbound cells with PBS and proceed to immunostaining (or continue the culture to obtain desired cell densities).


These cover-slips can also be used for staining suspension cells for microscopy, or for fixing suspension cells for confocal microscopy. What's more, you can also stimulate the cells while they are attached on the cover-slip, should you need to do so. You can work on non-adherent cells as if they were adherent!

Suspension cell immunofluorescence with Shi-fix™ coverslips


SIZE: ø 15mm x 0.17mm(#1.5)




Immunofluorescence analysis of Molt4 cells immobilized on SHI-FIX™ coverslips.



Immunofluorescence analysis of Daudi cells immobilized on SHI-FIX™ coverslips.



Immunofluorescence analysis of rat splenocytes immobilized on SHI-FIX™ coverslips



Immunofluorescence analysis of mouse thymocytes immobilized on SHI-FIX™ coverslips