



Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD, United Kingdom

everestbiotech.com

sales@everestbiotech.com

support@everestbiotech.com

Tel +44 1869 238326

Fax +44 1869 238327

**Research Use Only. Not for
diagnostic or therapeutic use.**

Storage: For long-term storage
keep aliquots at -20°C. (Store no
longer than 12 months at 4°C).
Minimize freezing and thawing.

EBP09223 - EBP09223 - Immunizing Peptide

Size: 100µg of dried peptide



Use

The peptide was used in the production of EB09223 – Goat Anti-ATG4D Antibody. It is currently untested in blocking applications.

Peptide Sequence

CSQLAQPTLRLPRT

Supplied As

100µg of dried peptide

Storage

Shipped at ambient temperature, store at -20°C.