



EB06718 - Goat Anti-Adiponectin Receptor 1 Antibody



Size: 100µg specific antibody in 200µl

This product is currently out of stock and is no longer available due to ongoing re-stocking problems. If re-stocking issues continue, Everest Biotech may have to take this product off the market very shortly. We apologize for any inconvenience caused. There are still some trial size samples available for the desperate customers.

Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD, United Kingdom

everestbiotech.com

sales@everestbiotech.com


support@everestbiotech.com

Tel +44 1869 238326

Fax +44 1869 238327

Research Use Only. Not for diagnostic or therapeutic use.

Storage: For long-term storage keep aliquots at -20°C. (Store no longer than 12 months at 4°C). Minimize freezing and thawing.



EB06718 (0.1 μ g/ml) staining of Human Muscle lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.