



Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD, United Kingdom

everestbiotech.com

sales@everestbiotech.com

support@everestbiotech.com

Tel +44 1869 238326

Fax +44 1869 238327

Research Use Only. Not for diagnostic or therapeutic use.

Storage: For long-term storage keep aliquots at -20°C. (Store no longer than 12 months at 4°C). Minimize freezing and thawing.

EB07240 - Goat Anti-AKR1C4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: type I 3-alpha-hydroxysteroid dehydrogenase, dihydrodiol dehydrogenase isozyme DD4, dihydrodiol dehydrogenase 4, chlordecone reductase, aldo-keto reductase family 1, member C4, 3-alpha-hydroxysteroid dehydrogenase variant, 3 alpha-hydroxysteroid dehydrogenase/dihydrodiol dehydrogenase 4, 3-alpha hydroxysteroid dehydrogenase, type I, MGC22581, HAKRA, DD4, CHDR, CDR, C11, 3-alpha-HSD, HGNC:387, aldo-keto reductase family 1, member C4, AKR1C4

Official Symbol: AKR1C4

Accession Number(s): NP_001809.2

Human GeneID(s): [1109](#)

Important Comments: This antibody may cross-react with AKR1C1

Immunogen

Peptide with sequence DPKYQRVELNDGH-C, from the N Terminus of the protein sequence according to NP_001809.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.


Western blot: Approx. 35kDa band observed in Human Liver lysates (calculated MW of 37.1kDa according to NP_001809.2). Recommended concentration: 0.1-0.3µg/ml.

IHC: In paraffin embedded Human Liver shows staining of distinct structures in the cytoplasm of hepatocytes. Recommended concentration, 2-4µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB07240 (0.1 μ g/ml) staining of human liver lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.



EB07240 (2.5 μ g/ml) staining of paraffin embedded Human Liver. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.