



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB11653 - Goat Anti-Antigen KI-67 / MKI67 (aa98-110) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** antigen identified by monoclonal antibody Ki-67, antigen KI-67, KIA, MKI67, OTTHUMP00000020739, OTTHUMP00000020740, proliferation-related Ki-67 antigen

**Official Symbol:** MKI67

**Accession Number(s):** NP\_002408.3; NP\_001139438.1

**Human GeneID(s):** [4288](#)

**Important Comments:** This antibody is expected to recognize both reported isoforms (NP\_002408.3; NP\_001139438.1).

### Immunogen

Peptide with sequence ENESLQNGRKSTE-C, from the N Terminus of the protein sequence according to NP\_002408.3; NP\_001139438.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:128000.

**Western blot:** Not yet tested - our routine western blotting protocol does not allow for the detection of proteins >250kDa. Therefore we cannot recommend an optimal concentration, and the antibody is an aspiring product.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human