



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12160 - Goat Anti-bcl-2-related protein A1 (rat) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Bcl2a1, BCL2-related protein A1, Bcl2a1d, B-cell leukemia/lymphoma 2 related protein A1, B-cell leukemia/lymphoma 2 related protein A1d, BCL2-related protein A1d, bcl-2-related protein A1

Official Symbol: Bcl2a1a

Accession Number(s): NP_596907.1

Non-Human GeneID(s): 12044 (mouse), 170929 (rat)

Important Comments: This antibody may cross-react with the almost identical proteins A1b (GeneID 12045) and A1d (GeneID 12047).

Immunogen

Peptide with sequence C-QKEVEKNLKPYLDD, from the internal region of the protein sequence according to NP_596907.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 19kDa band observed in Rat Spleen lysates (calculated MW of 20.4kDa according to NP_NP_596907.1). Recommended concentration: 0.1-0.3µg/ml. Antibody tested negative in mouse cells.

Species Reactivity

Tested: Rat

Expected from sequence similarity: Rat



EB12160 (0.1 μ g/ml) staining of Rat Spleen lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.