



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB12707 - Goat Anti-CNGB3 (aa175-186) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** CNGB3, cyclic nucleotide gated channel beta 3, ACHM1, CNG channel beta-3, cone photoreceptor cGMP-gated cation channel beta-subunit, cyclic nucleotide-gated cation channel beta-3, cyclic nucleotide-gated cation channel modulatory subunit

**Official Symbol:** CNGB3

**Accession Number(s):** NP\_061971.3

**Human GeneID(s):** [54714](#)

### Immunogen

Peptide with sequence C-ESDDKPTEHYR, from the internal region of the protein sequence according to NP\_061971.3.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested


**Peptide ELISA:** antibody detection limit dilution 1:8000.

**IHC:** In paraffin embedded Human Retina shows primarily staining of the pigmented epithelial cells. Recommended concentration: 2-4µg/ml.


### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human



EB12707 (2µg/ml) staining of paraffin embedded Human Retina. Microwaved antigen retrieval with citrate buffer pH 6, HRP-staining.



EB12707 Negative Control showing staining of paraffin embedded Human Retina, with no primary antibody.