



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12714 - Goat Anti-CYP2B6 (aa247-259) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CYP2B6, cytochrome P450, family 2, subfamily B, polypeptide 6, CPB6, CYP2B, CYP2B7, CYP2B7P, CYP2B6, EFVM, IIB1, P450, 1,4-cineole 2-exo-monooxygenase, cytochrome P450 2B6, cytochrome P450 IIB1, cytochrome P450, subfamily IIB (phenobarbital-inducible), polypeptide 6

Official Symbol: CYP2B6

Accession Number(s): NP_000758.1

Human GeneID(s): [1555](#)

Immunogen

Peptide with sequence C-HSVEKHRETLDPS, from the internal region of the protein sequence according to NP_000758.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 75kDa band observed in lysates of cell line HepG2 and K562 (calculated MW of 56.3kDa according to NP_000758.1). Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB12714 (0.1µg/ml) staining of HepG2 (A) and K562 (B) lysates (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.