



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB09921 - Goat Anti-DDA1 (aa32-42) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** C19orf58, cross-immune reaction antigen PCIA1, DET1 and DDB1 associated 1, MGC2594, PCIA1, DDA1

**Official Symbol:** DDA1

**Accession Number(s):** NP\_076955.1;

**Human GeneID(s):** [79016](#)

**Non-Human GeneID(s):** 66498 (mouse)

### Immunogen

Peptide with sequence C-PSVYLP TREYP, from the internal region of the protein sequence according to NP\_076955.1;

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Preliminary experiments gave no signal but low background in Human Umbilical Cord lysates and in lysates of cell line HEK293 at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Dog