



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09921 - Goat Anti-DDA1 (aa32-42) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: C19orf58, cross-immune reaction antigen PCIA1, DET1 and DDB1 associated 1, MGC2594, PCIA1, DDA1

Official Symbol: DDA1

Accession Number(s): NP_076955.1;

Human GeneID(s): [79016](#)

Non-Human GeneID(s): 66498 (mouse)

Immunogen

Peptide with sequence C-PSVYLP TREYP, from the internal region of the protein sequence according to NP_076955.1;

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave no signal but low background in Human Umbilical Cord lysates and in lysates of cell line HEK293 at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Dog