



Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD, United Kingdom

everestbiotech.com

sales@everestbiotech.com

support@everestbiotech.com

Tel +44 1869 238326

Fax +44 1869 238327

**Research Use Only. Not for
diagnostic or therapeutic use.**

Storage: For long-term storage
keep aliquots at -20°C. (Store no
longer than 12 months at 4°C).
Minimize freezing and thawing.

EB09141 - Goat Anti-DGCR8 / Pasha Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: pasha, Gy1, DGCRK6, C22orf12, DiGeorge syndrome critical region gene 8, DGCR8

Official Symbol: DGCR8

Accession Number(s): NP_073557.3

Human GeneID(s): [54487](#)

Non-Human GeneID(s): 94223 (mouse), 287954 (rat)

Immunogen

Peptide with sequence C-KRFDSEQVTVKKFRT, from the internal region of the protein sequence according to NP_073557.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave an approx. 75kDa band in Human Brain (amygdala) lysates after 0.2µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 86.0kDa according to NP_073557.3. The 75kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Cow