



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08152 - Goat Anti-DIAPH3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DIAPH3, diaphanous homolog 3 (Drosophila), RP11-26P21.1, DKFZP434C0931, DKFZp686A13178, DRF3, FLJ34705, diap3, diaphanous homolog 3, diaphanous-related formin 3

Official Symbol: DIAPH3

Accession Number(s): NP_001035982.1; NP_112194.2;
NP_001245295.1; NP_001245296.1; NP_001245297.1; NP_001245298.1

Human GeneID(s): [81624](#)

Immunogen

Peptide with sequence C-RTPMPKDVRQS, from the internal region (near the C terminus) of the protein sequence according to NP_001035982.1; NP_112194.2; NP_001245295.1; NP_001245296.1; NP_001245297.1; NP_001245298.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Species Reactivity

Tested:

Expected from sequence similarity: Human