



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06299 - Goat Anti-Dispatched homolog 1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DISPA, MGC13130, MGC16796, DKFZP434I0428, dispatched A, DISP1, dispatched homolog 1, dispatched homolog 1 (Drosophila), FLJ43740, MGC104180

Official Symbol: DISP1

Accession Number(s): NP_116279.2; NP_001337559.1

Human GeneID(s): [84976](#)

Immunogen

Peptide with sequence C-SELSGESLLIKTL, from the C Terminus of the protein sequence according to NP_116279.2; NP_001337559.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.


Immunofluorescence: Strong expression of the protein seen in the cell membrane of U87 cells. Recommended concentration: 2.5ug/ml.

Flow Cytometry: Flow cytometric analysis of HEK293 cells. Recommended concentration: 10ug/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Pig



EB06299 (2.5ug/ml) staining of U87. Detected by immunofluorescence. Data provided by Dr. Aashir Awan, University of Madrid, Spain.



EB06299 Flow cytometric analysis of paraformaldehyde fixed HEK293 cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control: Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.