



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08625 - Goat Anti-FPRL1 (aa13-26) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FPRL1, formyl peptide receptor-like 1, ALXR, FMLP-R-II, FMLPX, FPR2A, FPRH1, FPRH2, HM63, LXA4R, lipoxin A4 receptor (formyl peptide receptor related)

Official Symbol: FPRL1

Accession Number(s): NP_001453.1

Human GeneID(s): [2358](#)

Important Comments: Reported variants represent identical protein: NP_001005738.1 and NP_001453.1.

Immunogen

Peptide with sequence EEVSYESAGYTVLR-C, from the N Terminus of the protein sequence according to NP_001453.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary testing showed a band at approx 40kDa in Human Heart and approx. 37-38kDa in Human Spleen lysate after 1µg/ml antibody staining (calculated MW of 39.0kDa according to NP_001453.1). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Human