



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09584 - Goat Anti-FREM2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FREM2, FRAS1 related extracellular matrix protein 2, DKFZp686J0811, DKFZp781I048, KIAA1074, ECM3 homolog, FRAS1-related extracellular matrix protein 2, OTTHUMP00000018288, OTTHUMP00000018289

Official Symbol: FREM2

Accession Number(s): NP_997244.3

Human GeneID(s): [341640](#)

Immunogen

Peptide with sequence C-RLIRSEPTYNQPVQQ, from the internal region of the protein sequence according to NP_997244.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 351kDa according to NP_997244.3. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human