



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10727 - Goat Anti-Gsta3/Gsta5/Gst Yc2 (rat) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Gsta5, glutathione S-transferase Yc2 subunit, MGC112704, Yc2, GST A5-5, GST Yc2, glutathione S-transferase Yc-2, glutathione S-transferase alpha-5, glutathione S-transferase A3, Gsta3, glutathione S-transferase alpha 3, glutathione S-transferase alpha 3

Official Symbol: Gsta3

Accession Number(s): NP_001009920.1

Non-Human GeneID(s): 14859 (mouse), 494500 (rat)

Important Comments: Reported variants represent identical protein: NP_001009920.1, NP_001153211.1. Please note that this antibody is expected to cross-react to the homologous Mouse protein Gsta3

Immunogen

Peptide with sequence C-EMDPGIVDNFLLK, from the internal region of the protein sequence according to NP_001009920.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 25kDa band observed in Mouse and Rat lysates (calculated MW of 25.3kDa according to NP_001009920.1. Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Mouse, Rat

Expected from sequence similarity: Mouse, Rat



EB10727 (0.3 μ g/ml) staining of Rat Liver lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.