



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11744 - Goat Anti-Heme Oxygenase 1 (aa242-256) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: bK286B10, heat shock protein, 32-kD, heme oxygenase (decycling) 1, heme oxygenase 1, HMOX1, HO-1, HSP32, OTTHUMP00000198812

Official Symbol: HMOX1

Accession Number(s): NP_002124.1

Human GeneID(s): [3162](#)

Immunogen

Peptide with sequence C-NKVQDSAPVETPRGK, from the internal region of the protein sequence according to NP_002124.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 30kDa band observed in Human Liver lysates (calculated MW of 32.8kDa according to NP_002124.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB11744 (0.3 μ g/ml) staining of Human Liver lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.