



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09140 - Goat Anti-JAM2 / JAM-B Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: JAM2, junctional adhesion molecule 2, C21orf43, CD322, JAM-B, JAMB, PRO245, VE-JAM, VEJAM, JAM-IT/VE-JAM, OTTHUMP00000096100, junctional adhesion molecule B, vascular endothelial junction-associated molecule

Official Symbol: JAM2

Accession Number(s): NP_067042.1

Human GeneID(s): [58494](#)

Immunogen

Peptide with sequence C-NPRLGSQSTNSSYT , from the internal region of the protein sequence according to NP_067042.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave an approx 27kDa band in Human Heart, Liver and Placenta lysates after 0.2µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 33.2kDa according to NP_067042.1. The 27kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human