



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12057 - Goat Anti-LAT2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: LAT2, linker for activation of T cells family, member 2, LAB, NTAL, WBSCR15, WBSCR5, WSCR5, Williams-Beuren syndrome chromosomal region 15 protein, Williams-Beuren syndrome chromosomal region 5 protein, linker for activation of B-cells, linker for activation of T cells, transmembrane adaptor 2, linker for activation of T-cells family member 2, membrane-associated adapter molecule, non-T-cell activation linker

Official Symbol: LAT2

Accession Number(s): NP_054865.2

Human GeneID(s): [7462](#)

Important Comments: Reported variants represent identical protein: NP_115853.2, NP_115852.1, NP_054865.2

Immunogen

Peptide with sequence C-RSEKIYQQRSLRED, from the internal region of the protein sequence according to NP_054865.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx. 26kDa band observed in Human Peripheral Blood Monocytes lysates and Human Spleen, and approx. 30kDa in lysates of cell lines Daudi and U937 (calculated MW of 26.6kDa according to NP_054865.2). Recommended concentration: 0.1-2µg/ml. Primary incubation 1 hour at room temperature.

Immunofluorescence: Strong expression of the protein seen in U937 cells.


Recommended concentration: 10µg/ml.

Flow Cytometry: Flow cytometric analysis of U937 cells. Recommended concentration: 10ug/ml.


Species Reactivity

Tested: Human


Expected from sequence similarity: Human




EB12057 (0.1 μ g/ml) staining of Daudi (A), U937 (B) and (2 μ g/ml) PBM (C) cell lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB12057 (0.3 μ g/ml) staining of Human Spleen lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB12057 Immunofluorescence analysis of paraformaldehyde fixed U937 cells immobilized on ShifixTM coverslip, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing strong membrane and some cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).



EB12057 Flow cytometric analysis of paraformaldehyde fixed U937 cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control: Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.