



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07794 - Goat Anti-LNX1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: LNX1, ligand of numb-protein X 1, LNX, MPDZ, PDZRN2, multi-PDZ-domain-containing protein

Official Symbol: LNX1

Accession Number(s): NP_116011.1; NP_001119800.1

Human GeneID(s): [84708](#)

Immunogen

Peptide with sequence C-PDAYRPRDDSFH, from the internal region of the protein sequence according to NP_116011.1; NP_001119800.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave next to 70kDa bands at approx 75kDa and 27kDa in Human Kidney and Lung lysates after 0.5 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the extra bands we observe given the calculated size of 69.7kDa according to NP_116011.1. The detected bands were successfully blocked by incubation with the immunizing peptide (and BLAST results with the immunizing peptide sequence did not identify any other proteins to explain the additional bands). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human