



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08496 - Goat Anti-Lubricin / PRG4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PRG4, proteoglycan 4, CACP, FLJ32635, HAPO, JCAP, MSF, SZP, bG174L6.2, Jacobs camptodactyly-arthropathy-pericarditis syndrome, articular superficial zone protein, bG174L6.2 (MSF: megakaryocyte stimulating factor), camptodactyly, arthropathy, coxa vara, pericarditis syndrome, lubricin, megakaryocyte stimulating factor, proteoglycan 4, (megakaryocyte stimulating factor, articular superficial zone protein, camptodactyly, arthropathy, coxa vara, pericarditis syndrome)

Official Symbol: PRG4

Accession Number(s): NP_005798.3; NP_001121180.2; NP_001290161.1

Human GeneID(s): [10216](#)

Immunogen

Peptide with sequence DYESFCAEVHNP, from the internal region of the protein sequence according to NP_005798.3; NP_001121180.2; NP_001290161.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Species Reactivity

Tested:

Expected from sequence similarity: Human