



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09331 - Goat Anti-MRGPRX3 / MRGX3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MRGPRX3, MAS-related GPR, member X3, GPCR, MRGX3, SNSR1, G protein-coupled receptor MRGX3, G protein-coupled receptor SNSR1, G protein-coupled receptor SNSR2, mas-related GPCR member X3

Official Symbol: MRGPRX3

Accession Number(s): NP_473372.2

Human GeneID(s): [117195](#)

Important Comments: This antibody may cross-react with human MRGPRX1 and X4 (GeneID: 259249, 117196; MW: 36.25kDa and 36.43kDa).

Immunogen

Peptide with sequence C-QLPQETLELSGSRL, from the Internal region of the protein sequence according to NP_473372.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product.

Species Reactivity

Tested:

Expected from sequence similarity: Human