



Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD, United Kingdom

everestbiotech.com

sales@everestbiotech.com

support@everestbiotech.com

Tel +44 1869 238326

Fax +44 1869 238327

Research Use Only. Not for diagnostic or therapeutic use.

Storage: For long-term storage keep aliquots at -20°C. (Store no longer than 12 months at 4°C). Minimize freezing and thawing.

EB05238 - Goat Anti-MTMR1 (isoform 1) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: myotubularin-related protein 1, myotubularin related protein 1, MTMR1

Official Symbol: MTMR1

Accession Number(s): NP_003819.1

Human GeneID(s): [8776](#)

Important Comments: Please note this will only recognise Human isoform 1 (NP_003819.1), not isoform 2 (NP_789746).

Immunogen

Peptide with sequence C-SSPSHSATS VHTSV, from the C Terminus of the protein sequence according to NP_003819.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave no signal but low background in Human Liver and Placenta lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human