

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11103 - Goat Anti-Muc16 (Rat) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: 1110008I14Rik, CA125, CA-125, CA125 ovarian cancer antigen, Gm1480, MUC-16, mucin 16, cell surface associated, mucin-16, ovarian cancer-related tumor marker CA125, ovarian carcinoma antigen CA125, RGD1308422, MUC16

Official Symbol: Muc16

Accession Number(s): XP_235886.5

Non-Human GeneID(s): 73732 (mouse), 315451 (rat)

Immunogen

Peptide with sequence C-QRHRLAYLSHLDLRK, from the internal region of the protein sequence according to XP_235886.5.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 794kDa according to XP_235886.5. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Rat, Mouse