



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11673 - Goat Anti-Nbeal2 (mouse aa543-556) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: 1110014F23Rik, BC042396, mKIAA0540, Nbeal2, neurobeachin-like 2, neurobeachin-like protein 2, OTTMUSP00000043650, OTTMUSP00000043652

Official Symbol: Nbeal2

Accession Number(s): NP_899099.2

Non-Human GeneID(s): 235627 (mouse)

Immunogen

Peptide with sequence C-DSESSGNESQKARH, from the internal region of the protein sequence according to NP_899099.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 303kDa according to NP_899099.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Mouse