



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB06284 - Goat Anti-P8 / NUPR1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** NUPR1, Nuclear Protein 1, P8, COM1, p8 protein (candidate of metastasis 1), p8 protein, nuclear protein, transcriptional regulator, 1, candidate of metastasis 1

**Official Symbol:** NUPR1

**Accession Number(s):** NP\_001035948.1; NP\_036517

**Human GeneID(s):** [26471](#)

### Immunogen

Peptide with sequence ATFPATSAPQQP-C, from the N Terminus of the protein sequence according to NP\_001035948.1; NP\_036517.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Preliminary experiments gave no signal but low background in HepG2 and Human Pancreas extracts at up to 1µg/ml. Please note that we would not expect our routine western blot protocol to reliably detect proteins as small as the predicted size of this protein (approx 9kDa according to NP\_036517). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human