



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12163 - Goat Anti-Prdm9 (mouse, aa440-452) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Prdm9, PR domain containing 9, RP23-64O5.3, BC012016, Dsbc1, G1-419-29, Meisetz, Rcr1, repro7, PR domain zinc finger protein 9, PR domain-containing protein 9, PR-domain zinc finger protein, VPR domain containing 9, histone-lysine N-methyltransferase PRDM9, hybrid sterility protein 1, meiosis-induced factor containing a PR/SET domain and zinc-finger motif

Official Symbol: Prdm9

Accession Number(s): NP_659058.2

Non-Human GeneID(s): 213389 (mouse)

Immunogen

Peptide with sequence C-QEQHVDSQNKNDK, from the internal region of the protein sequence according to NP_659058.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Preliminary experiments gave an approx 70kDa band in Mouse fetal Brain lysates and in lysates of cell line NIH3T3 after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 97.8kDa according to NP_659058.2. The 70kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Mouse