



Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD, United Kingdom

everestbiotech.com

sales@everestbiotech.com

support@everestbiotech.com

Tel +44 1869 238326

Fax +44 1869 238327

Research Use Only. Not for diagnostic or therapeutic use.

Storage: For long-term storage keep aliquots at -20°C. (Store no longer than 12 months at 4°C). Minimize freezing and thawing.

EB06566 - Goat Anti-Prostaglandin dehydrogenase 1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: hydroxyprostaglandin dehydrogenase 15-(NAD), 15-PGDH, PGDH1, HPGD, Prostaglandin dehydrogenase 1

Official Symbol: HPGD

Accession Number(s): NP_000851.2

Human GeneID(s): [3248](#)

Immunogen

Peptide with sequence C-DYDTTPFQAKTQ, from the C Terminus of the protein sequence according to NP_000851.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.


Western blot: Experiments gave bands at approx 26kDa and 50kDa in Human Duodenum lysates after 0.1µg/ml antibody staining. These bands correspond to earlier findings in literature with different antibodies (PMID: 11889207). This protein has a calculated MW of 29.0kDa according to NP_000851. Recommended concentration: 0.1-1µg/ml.

IHC: In paraffin embedded Human Colon shows strong cytoplasm staining of enterocytes. Recommended concentration, 2-4µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB06566 (0.1µg/ml) staining of Human Duodenum lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.



EB06566 (2.5µg/ml) staining of paraffin embedded Human Colon. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.