

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08277 - Goat Anti-PSD3 (aa 395-409) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PSD3, pleckstrin and Sec7 domain containing 3, DKFZp761K1423, EFA6R, HCA67, ADP-ribosylation factor guanine nucleotide factor 6, hepatocellular carcinoma-associated antigen 67

Official Symbol: PSD3

Accession Number(s): NP_056125.3; NP_996792.1

Human GeneID(s): [23362](#)

Important Comments: This antibody is expected to recognise the reported isoform (NP_056125.3).

Immunogen

Peptide with sequence C-QENKQHLEKTPKPER, from the internal region of the protein sequence according to NP_056125.3; NP_996792.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Preliminary experiments in Human Brain (Cerebellum and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).

Species Reactivity

Tested:

Expected from sequence similarity: Human