



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10368 - Goat Anti-PSMB9 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: beta1i, LMP2, low molecular mass protein 2, macropain chain 7, MGC70470, multicatalytic endopeptidase complex chain 7, proteasome (prosome, macropain) subunit, beta type, 9 (large multifunctional peptidase 2), proteasome beta 9 subu, PSMB6i, RING12, PSMB9

Official Symbol: PSMB9

Accession Number(s): NP_002791.1

Human GeneID(s): [5698](#)

Non-Human GeneID(s): 16912 (mouse), 24967 (rat)

Immunogen

Peptide with sequence C-RNISYKYRED, from the internal region of the protein sequence according to NP_002791.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:2000.

Western blot: Approx 22kDa band observed in Human Thymus lysates and approx. 20kDa in Rat Thymus (calculated MW of 23.3kDa according to Human NP_002791.1 and Rat NP_036840.2). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow, Pig



EB10368 (1 μ g/ml) staining of Human Thymus lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.



EB10368 (0.3 μ g/ml) staining of Rat Thymus lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.