



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07411 - Goat Anti-Quaking / QKI Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: QKI, quaking homolog, KH domain RNA binding, DKFZp586I0923, Hqk, QK, QK3, RNA binding protein HQK, homolog of mouse quaking QKI (KH domain RNA binding protein), quaking homolog, KH domain RNA binding

Official Symbol: QKI

Accession Number(s): NP_006766.1; NP_996735.1; NP_996736.1; NP_996737.1

Human GeneID(s): [9444](#)

Non-Human GeneID(s): 19317 (mouse)

Important Comments: This antibody is expected to recognise all reported isoforms (NP_006766.1; NP_996735.1; NP_996736.1; NP_996737.1).

Immunogen

Peptide with sequence ETKEKPKPTPDY-C, from the N Terminus of the protein sequence according to NP_006766.1; NP_996735.1; NP_996736.1; NP_996737.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Approx 38kDa band observed in mouse brain lysates (calculated MW of 35.8kDa according to NP_068681.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB07411 (0.3 μ g/ml) staining of mouse brain lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.