



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06243 - Goat Anti-RANBP16 / Exportin 7 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: XPO7, Exportin 7, RANBP16, KIAA0745, RAN binding protein 16

Official Symbol: XPO7

Accession Number(s): NP_001093631.1; NP_055839; NP_001093632.1

Human GeneID(s): [23039](#)

Non-Human GeneID(s): 65246 (mouse)

Important Comments: This antibody is expected to recognise all the reported isoforms (NP_001093631.1; NP_055839; NP_001093632.1).

Immunogen

Peptide with sequence C-NSTYGVSNDMMMS, from the C Terminus of the protein sequence according to NP_001093631.1; NP_055839; NP_001093632.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 110kDa band observed in Hela lysates (predicted size of approx. 124kDa according to NP_055839). Recommended concentration: 0.5-1µg/ml

IHC: Paraffin embedded Human Testis. Recommended concentration: 5µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Dog, Pig, Cow



EB06243 staining (0.5µg/ml) of Hela lysate (RIPA buffer, 35µg total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.



EB06243 (5µg/ml) staining of paraffin embedded Human Testis. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.