

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07208-T - Goat Anti-SOD1 Antibody - Trial

Size: 20µg specific antibody in 40µl



Target Protein

Principal Names: SOD1, superoxide dismutase 1, soluble (amyotrophic lateral sclerosis 1 (adult)), HGNC:11179, ALS, ALS1, IPOA, SOD, Cu/Zn superoxide dismutase, SOD, soluble, indophenoloxidase A, superoxide dismutase 1, soluble, superoxide dismutase, cystolic, HEL-S-44, hSod1, epididymis secretory protein Li 44

Official Symbol: SOD1

Accession Number(s): NP_000445.1

Human GeneID(s): [6647](#)

Non-Human GeneID(s): 20655 (mouse), 24786 (rat)

Immunogen

Peptide with sequence C-SRKHGGPKDEERH, from the internal region of the protein sequence according to NP_000445.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 18kDa band observed in Mouse Brain and Rat Spinal Cord lysates, and in lysates of cell line NIH3T3. Approx 20kDa observed in lysates of cell lines HEK293, HepG2 and MCF7 (calculated MW of 15.9kDa according to Human NP_000445.1, Mouse NP_035564.1 and Rat NP_058746.1). Recommended concentration: 0.01-0.03µg/ml. Primary incubation 1 hour at room temperature. This antibody has been successfully used in Western blot on Mouse: Frazziano G et al. (2014) PMID: 24213612.

Immunofluorescence: Strong expression of the protein seen in the nucleus and cytoplasm of U2OS cells and the cytoplasm of A431 cells. Recommended concentration: 10µg/ml.

Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human, Mouse, Rat, Dog

Specific References

This antibody has been successfully used in the following paper:

Krzysztof Sikorski, Adi Mehta, Marit Inngjerdningen, Flourina Thakor, Simon Kling, Tomas Kalina, Tuula A. Nyman, Maria Ekman Stensland, Wei Zhou, Gustavo A. De Souza, Lars Holden, Jan Stuchly, Markus Templin and Fridtjof Lund-Johansen
A high-throughput pipeline for validation of antibodies
Nat Methods. 2018 Nov;15(11):909-912
PMID: 30377371


This antibody has been successfully used in Western blot on Mouse:

Frazziano G, Al Ghoulé I, Baust J, Shiva S, Champion HC, Pagano PJ.


Nox-derived ROS are acutely activated in pressure overload pulmonary hypertension: indications for a seminal role for mitochondrial Nox4.

Am J Physiol Heart Circ Physiol. 2014 Jan 15;306(2):H197-205.


PMID: 24213612




EB07208 (0.02 μ g/ml) staining of Mouse Brain lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB07208 (0.01 μ g/ml) staining of NIH3T3 (A), HEK293 (B), HepG2 (C) and MCF7 (D) lysates (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB07208 Immunofluorescence analysis of paraformaldehyde fixed U2OS cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing nuclear and cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).



EB07208 Immunofluorescence analysis of paraformaldehyde fixed A431 cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).