



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06180-T - Goat Anti-Syntrophin (alpha 1) Antibody - Trial

Size: 20µg specific antibody in 40µl



Target Protein

Principal Names: SNTA1, syntrophin, alpha 1 (dystrophin-associated protein A1, 59kDa, acidic component), SNT1, TACIP1, dJ1187J4.5, pro-TGF-alpha cytoplasmic domain-interacting protein 1, dystrophin-associated protein A1, 59kD, acidic component, syntrophin, alpha 1 (dystrophin-associated protein A1, 59kD, acidic component), LQT12, acidic alpha 1 syntrophin

Official Symbol: SNTA1

Accession Number(s): NP_003089.1

Human GeneID(s): [6640](#)

Non-Human GeneID(s): 20648 (mouse)

Immunogen

Peptide with sequence ASGRRAPRTGLLE-C, from the N Terminus of the protein sequence according to NP_003089.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 60kDa band observed in Human and Mouse Muscle lysates and in rodent Heart lysates (calculated MW of 53.9kDa according to Human NP_003089.1).

Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human, Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB06180 (0.01 µg/ml) staining of Human Muscle lysate (35 µg protein in RIPA buffer) with (B) and without (A) blocking with the immunising peptide. Primary incubation was 1 hour. Detected by chemiluminescence.