



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05984 - Goat Anti-USP11 / UHX1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: USP11, UHX1, ubiquitin specific protease 11, deubiquitinating enzyme 11, ubiquitin thiolesterase 11, ubiquitin carboxyl-terminal hydrolase, X-linked, ubiquitin specific peptidase 11, RP4-659F15.2, ubiquitin-specific processing protease 11

Official Symbol: USP11

Accession Number(s): NP_004642.2

Human GeneID(s): [8237](#)

Immunogen

Peptide with sequence CSSPPSSEFMDVN, from the C Terminus of the protein sequence according to NP_004642.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 120kDa band observed in lysates of cell line Jurkat (calculated MW of 110kDa according to NP_004642.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Ideguchi et al, Biochem J. 2002 Oct 1;367(Pt 1):87-95.; PMID: 12084015). Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB05984 staining (1 μ g/ml) of Jurkat lysate (RIPA buffer, 30 μ g total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.