



Macrophage M1 + M2

Macrophages are a highly diverse cell population, exhibiting phenotypic heterogeneity and functional specialization, and play crucial roles in both innate and adaptive immunity. There are two main subtypes: M1 (classically activated, exhibiting a pro-inflammatory phenotype and are necessary for pathogenic defence) and M2 (alternatively activated, exhibiting an anti-inflammatory phenotype, involved in the resolution of inflammation and tissue repair).


Please view our M1 + M2 Marker antibodies below:



EB12699 (5µg/ml) staining of paraffin embedded Human Spleen. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



EB12715 Immunofluorescence analysis of paraformaldehyde fixed THP-1 cells immobilized on Shi-fix™ plus cover-slips.



EB07827 Flow cytometric analysis of A549 cells. Recommended concentration: 10µg/.ml.

Everest code	Principal name	Applications Tested	Species reactivity	Size available
EB07691	Goat Anti-AKT1 Antibody	Pep-ELISA, WB, IF, FC	Human, Mouse	100µg/200µl
EB06875	Goat Anti-AKT1 Antibody	Pep-ELISA, WB	Human, Mouse, Cow	100µg/200µl
EB06772	Goat Anti-AKT2 Antibody	Pep-ELISA, WB, IHC	Human, Mouse, Rat, Dog	100µg/200µl
EB07159	Goat Anti-AKT2 Antibody	Pep-ELISA, IF, FC	Human, Mouse, Rat	100µg/200µl
EB06562	Goat Anti-AKT3 Antibody	Pep-ELISA, WB, IHC, IF, FC	Human, Mouse, Rat, Dog, Cow	100µg/200µl
EB06795	Goat Anti-Arginase I Antibody	Pep-ELISA, WB, IF	Human, Dog, Cow, Pig	100µg/200µl
EB07706	Goat Anti-Arginase, type 1 / ARG1 Antibody	Pep-ELISA, WB, IHC, IF	Human, Dog	100µg/200µl
EB07707	Goat Anti-Arginase, type 1 / arg1(rat) Antibody	Pep-ELISA, WB, IF	Mouse, Rat	100µg/200µl
EB09667	Goat Anti-CCL20 Antibody	Pep-ELISA, IHC	Human, Cow, Pig	100µg/200µl
EB09000	Goat Anti-CCL3L1 Antibody	Pep-ELISA	Human	100µg/200µl

Everest code	Principal name	Applications Tested	Species reactivity	Size available
EB07239	Goat Anti-CCR2 Antibody	Pep-ELISA, IF, IHC, FC	Human	100µg/200µl
EB12341	Goat Anti-CD163 (aa871-882) Antibody	Pep-ELISA, WB, IHC	Human	100µg/200µl
EB12340	Goat Anti-CD163 Anti-body	Pep-ELISA, WB, IHC, FC	Human	100µg/200µl
EB06621	Goat Anti-CD274 / PD-L1 Antibody	Pep-ELISA, WB, IF, FC	Human	100µg/200µl
EB12699	Goat Anti-CD68 Antibody	Pep-ELISA, WB, IHC, FC	Human	100µg/200µl
EB09701	Goat Anti-CD80 Antibody	Pep-ELISA, WB, IHC, FC	Human	100µg/200µl
EB10738	Goat Anti-CD86 / B7-2 (aa273-286) Antibody	Pep-ELISA	Human	100µg/200µl
EB10739	Goat Anti-CD86 / B7-2 (aa73-87) Antibody	Pep-ELISA, IHC	Human	100µg/200µl
EB10257	Goat Anti-CSF1R (aa 700-710) Antibody	Pep-ELISA, WB-Trf	Human, Mouse, Rat, Dog, Cow	100µg/200µl
EB10258	Goat Anti-CSF1R (aa734-744) Antibody	Pep-ELISA, FC, IHC, WB	Human, Dog, Cow	100µg/200µl
EB12715	Goat Anti-DC-SIGN / CD209 Antibody	Pep-ELISA, WB, IHC, IF, FC	Human	100µg/200µl
EB12555	Goat Anti-FCGR1A/CD64 Antibody	Pep-ELISA, WB, FC	Human	100µg/200µl
EB12342	Goat Anti-IL-1 beta Anti-body	Pep-ELISA, WB	Human	100µg/200µl
EB10107	Goat Anti-IL1R2 Antibody	Pep-ELISA, WB, IF, FC	Human	100µg/200µl
EB11755	Goat Anti-IL-6 (aa138-148) Antibody	Pep-ELISA, IF, IHC	Human	100µg/200µl
EB11331	Goat Anti-IRF4 (aa393-407) Antibody	Pep-ELISA, WB, IHC, FC	Human, Mouse, Rat, Dog, Cow	100µg/200µl
EB11330	Goat Anti-IRF4 Antibody	Pep-ELISA, WB	Human	100µg/200µl
EB07825	Goat Anti-ITGAM / CD11B Antibody	Pep-ELISA, IF, IHC	Human	100µg/200µl
EB07241	Goat Anti-PPARG Anti-body	Pep-ELISA, IF	Human, Mouse, Rat	100µg/200µl
EB11961	Goat Anti-RANTES Anti-body	Pep-ELISA, IHC	Human	100µg/200µl
EB05040	Goat Anti-SOCS1 Antibody	Pep-ELISA, FC, IF, IHC	Human, Mouse, Dog, Cow, Rat, Pig	100µg/200µl
EB05104	Goat Anti-SOCS2 Antibody	Pep-ELISA, IHC	Human, Mouse, Rat	100µg/200µl
EB07852	Goat Anti-TGM2 Antibody	Pep-ELISA, FC, WB	Human, Dog, Cow	100µg/200µl
EB07827	Goat Anti-TGM2 Antibody	Pep-ELISA, WB, IHC, FC	Human	100µg/200µl
EB11080	Goat Anti-Transglutaminase 2 (aa17-29) Antibody	Pep-ELISA, IF	Human, Mouse, Rat, Pig, Cow	100µg/200µl
EB11081	Goat Anti-Transglutaminase 2 (aa578-591) Antibody	Pep-ELISA, WB	Human, Mouse, Rat, Pig, Cow	100µg/200µl