


GOAT ANTI-LAMIN A / LAMIN A/C ANTIBODY

SKU: EB06696



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names Charcot-Marie-Tooth disease, axonal, type 2B1|progeria 1 (Hutchinson-Gilford type)|70 kDa lamin|lamin A/C|lamin A/C|lamin A|LGMD1B|CMT2B1|CMD1A|PRO1|LMNC|LMN1|LDP1|HGPS|FPLD|EMD2|LFP|FPL|LMNA

Accession ID NP_733821.1; NP_733822.1

Blocking Peptide EBP06696

Immunogen Peptide with sequence PRTQSPQNCSIM, from the C Terminus of the protein sequence according to NP_733821.1; NP_733822.1.

Product Comments This antibody is expected to recognise isoforms 1 and 3 (NP_733821.1; NP_733822.1).

Peptide Sequence PRTQSPQNCSIM

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	4000
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Prostate. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

GALLERY IMAGES

