


GOAT ANTI-RSK2 ANTIBODY

SKU: EB06985



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias ribosomal protein S6 kinase, 90kD, polypeptide 3|mental retardation, X-linked 19|insulin-stimulated protein kinase 1|RSK|CLS|RP11-393H10.3|pp90RSK2|p90-RSK2|S6K-alpha3|RSK2|MRX19|MAPKAPK1B|ISPK-1|HU-3|HGNC:10432|ribosomal protein S6 kinase, 90kDa, polypeptide 3|RPS6KA3

Accession ID NP_004577.1

Immunogen Peptide with sequence C-DDESQAMQTVGVH, from the internal region of the protein sequence according to NP_004577.1.

Peptide Sequence C-DDESQAMQTVGVH

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat

Reactive Species Human

Human Gene ID 6197

Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Skin. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, IHC

GALLERY IMAGES

