

EBP-PAIR-001-HRP - IMMUNIZING PEPTIDE

SKU: EBP-PAIR-001-HRP



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	P0DTC2; NC_045512.2/MN908947; P59594
Immunogen	Peptide with sequence PROPRIETARY'-C, from the of the protein sequence according to P0DTC2; NC_045512.2/MN908947; P59594.
Product Comments	The peptide was used in the production of EB-PAIR-001-HRP - Anti-SARS-CoV2 Spike protein S1 Matched ELISA pair Antibody. It is currently untested in blocking applications.
Peptide Sequence	PROPRIETARY'-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png