

EBP05007 - IMMUNIZING PEPTIDE

SKU: EBP05007



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_006143.2 |
| Immunogen | Peptide with sequence C-PFTRLRHNGPPPV, from the C Terminus of the protein sequence according to NP_006143.2. |
| Product Comments | The peptide was used in the production of EB05007 - Goat Anti-LRMP / JAW1 Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | C-PFTRLRHNGPPPV |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |