

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP05010 - IMMUNIZING PEPTIDE**

**SKU:** EBP05010

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 003870.4; NP 001189445.1; NP 001189446.1; NP 001189444.1 ID

Peptide with sequence C-QHTLRKKLILSYT, from the C Terminus of the protein sequence according to Immunogen

NP 003870.4; NP 001189445.1; NP 001189446.1; NP 001189444.1.

The peptide was used in the production of EB05010 - Goat Anti-CFLAR / FLIP Antibody. It is currently untested **Product** 

**Comments** in blocking applications.

**Peptide** 

C-QHTLRKKLILSYT Sequence

Shipping

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png