

EBP05010 - IMMUNIZING PEPTIDE

SKU: EBP05010



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003870.4; NP_001189445.1; NP_001189446.1; NP_001189444.1
Immunogen	Peptide with sequence C-QHTLRKKLILSYT, from the C Terminus of the protein sequence according to NP_003870.4; NP_001189445.1; NP_001189446.1; NP_001189444.1.
Product	The peptide was used in the production of EBP05010 - Goat Anti-CFLAR / FLIP Antibody. It is currently untested
Comments	in blocking applications.
Peptide Sequence	C-QHTLRKKLILSYT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png