

EBP05017 - IMMUNIZING PEPTIDE

SKU: EBP05017



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002036
Immunogen	Peptide with sequence C-GKEEPEADQEHA, from the C Terminus of the protein sequence according to NP_002036.
Product	The peptide was used in the production of EB05017 - Goat Anti-GAP43 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-GKEEPEADQEHA
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png