

EBP05028 - IMMUNIZING PEPTIDE

SKU: EBP05028



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002620.1; NP_000281.2; NP_001025062.1
Immunogen	Peptide with sequence C-KAMEAVAAQGKAKK, from the C Terminus of the protein sequence according to NP_002620.1; NP_000281.2; NP_001025062.1.
Product	The peptide was used in the production of EB05028 - Goat Anti-PGAM1 / PGAM2 / PGAM4 Antibody. It is
Comments	currently untested in blocking applications.
Peptide Sequence	C-KAMEAVAAQGKAKK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)