

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05030 - IMMUNIZING PEPTIDE

SKU: EBP05030

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 002856.1 ID

Peptide with sequence HAGNQGGQQAGGGCC, from the C Terminus of the protein sequence according to **Immunogen**

NP_002856.1.

The peptide was used in the production of EB05030 - Goat Anti-RAB2 / RAB2A Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

HAGNQGGQQAGGGCC Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png