

EBP05036 - IMMUNIZING PEPTIDE

SKU: EBP05036



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001008707.1; NP_004425.2;
Immunogen	Peptide with sequence C-GGKDTSIMQWRVI, from the C Terminus of the protein sequence according to NP_001008707.1; NP_004425.2;.
Product Comments	The peptide was used in the production of EBP05036 - Goat Anti-EML1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GGKDTSIMQWRVI
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiling_medium.png