

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05037 - IMMUNIZING PEPTIDE

SKU: EBP05037

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_036287.1; NP_001180197.1; NP_001180198.1 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \begin{array}{l} \text{Peptide with sequence C-GKDTSVLQWRVV, from the C Terminus of the protein sequence according to} \\ \text{NP}_036287.1; \ \text{NP}_001180197.1; \ \text{NP}_001180198.1. \end{array}$

The peptide was used in the production of EB05037 - Goat Anti-EML2 / ELP70 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-GKDTSVLQWRVV Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade