

EBP05037 - IMMUNIZING PEPTIDE

SKU: EBP05037



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_036287.1; NP_001180197.1; NP_001180198.1
Immunogen	Peptide with sequence C-GKDTSVLQWRVV, from the C Terminus of the protein sequence according to NP_036287.1; NP_001180197.1; NP_001180198.1.
Product Comments	The peptide was used in the production of EBP05037 - Goat Anti-EML2 / ELP70 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GKDTSVLQWRVV
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)