



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05040 - IMMUNIZING PEPTIDE

SKU: EBP05040

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_003736.1 ID

Immunogen Peptide with sequence C-VLRDYLSSFPFQI, from the C Terminus of the protein sequence according to NP_003736.1.

Product The peptide was used in the production of EB05040 - Goat Anti-SOCS1 Antibody. It is currently untested in blocking

Comments applications.

Peptide

C-VLRDYLSSFPFQI Sequence

Shipping Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade